

# Emergency-Ventilation with the OXYLATOR®

## OXYLATOR® EMX



- Open the main valve of the cylinder.
- Rotate the airway-pressure selector to 20 cmH<sub>2</sub>O.



- Press the O<sub>2</sub> button and turn it clockwise
- Now the OXYLATOR® is in automatic mode



- Place the mask on the patient's face.
- Assure that there is a proper mask seal.
- Hold the mask with both hands and tilt the head.



- Manual mask ventilation – press the O<sub>2</sub> button for inspiratory flow, release it for passive exhalation.
- Every second of inspiratory time delivers 500 ml of volume of oxygen towards the patient's airway.
- Try to achieve an inspiratory time of 1.5 to 2 seconds not exceeding the pressure setting of 20 cmH<sub>2</sub>O.



- If the patient is intubated, you can adjust the pressure setting between 20 – 45 cmH<sub>2</sub>O to achieve 1.5 to 2 seconds of inspiratory time.



- Manual ventilation with an intubated patient - press the O<sub>2</sub> button for inspiratory flow, release it for passive exhalation.
- Every second of inspiratory time delivers 500 ml of volume of oxygen towards the patient's airway.



The addition of a capnometer or an ETCO<sub>2</sub> monitoring device to the OXYLATOR® allows you to adjust the OXYLATOR® to the patient's needs.